

VISION™

Eye Wash Station with 500ml Eye Cap Bottles



**Robust eye wash station.
Wipe clean, hardwearing, ABS plastic.
Flourescent, printed instructions.
Includes integral eye bath.**

Dual purpose HSE compliant eye wash station. Wipeable, hard-wearing ABS plastic with flourescent printed instructions; useful if a power cut happens so the station can be located. Dust proof heavy duty cover keeps contents clean. Central Mirror for easy application. Contains easy-to-squeeze bottle with an eye cap, allowing the eye to be irrigated easily and efficiently.



**Eye wash bottle with
integral eye bath.**

Vision Eye Wash Station Green/Clear Lid

REF	PRODUCT	SIZE	UNIT	PACK
919	Vision Eye Wash Station	36cmH x 51cmW x 10cmD	1	1
918	Vision Eye Wash Station - empty	36cmH x 51cmW x 10cmD	1	1
933	Refill for Vision Eye Wash Station	-	1	1

CONTENTS

2	Dressings, Eye Pad	5	Eye Wash, Pods 20ml
2	Eye Wash, Bottles 500ml	1	Mirror

